TransBioLine

Translational Safety Biomarker Pipeline

821283 – TransBioLine

Translational Safety Biomarker Pipeline

WP7 – Assay and Sample Management

D7.02 SOP for sample collection, processing, storage, and shipment

Lead contributor	Michael Hummel (15 – Charité), Dana Briesemeister (15 – Charité)
Other contributors	Organ WPs, 9- TamiRNA, 10 – SIGNATOPE, Humaira Mahmood (21 - PFIZER), Jane Grove (6 – UNOTT)

Due date	31.05.2019
Delivery date	22.01.2020
Deliverable type	R
Dissemination level	PU

Description of Action	Version	Date
	V1	11/09/2017

Acknowledgement

This project has received funding from the Innovative Medicines Initiative 2 Joint Undertaking under grant agreement No 821283. This Joint Undertaking receives support from the European Union's Horizon 2020 research and innovation programme and EFPIA companies.

Index

Document History	2
Definitions	3
Publishable Summary	4
1. Introduction	5
2. Methods	5
2.1 Retrieval of information	5
2.1.1 Request of information for analyses	5
2.1.2 Information study protocols/SOPs	5
2.2. Harmonization of SOPs	5
2.2.1 SOP for sample collection and processing	5
3. Results	6
3.1 SOP for sampling, processing and intermediate storage	6
3.2. SOP for storage	6
3.2.3 SOP for shipment	6
4. Discussion	7
5. Conclusion	7
Annexes	7
Annex 1	7

Document History

Version	Date	Description
V1.0	20.05.2019	First Draft
V2.0	05.11.2019	Second Draft
V3.0	31.12.2019	Comments
V4.0	06.01.2020	Draft
V5.0	22.01.2020	Final Version

Definitions

Participants of the TransBioLine Consortium are referred to herein according to the following codes:

UZH. University of Zurich. Switzerland

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UNOTT. The University of Nottingham. United Kingdom

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MetaHeps. MetaHeps GmbH. Germany

APHP. Assistance Publique - Hopitaux de Paris. France

SYNAPSE. Synapse Research Management Partners S.L. Spain

Charité. Charité - Universitätsmedizin Berlin. Germany

UNEW. University of Newcastle Upon Tyne. United Kingdom

ULIV. The University of Liverpool. United Kingdom

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Grant Agreement. The agreement signed between the beneficiaries and the IMI JU for the undertaking of the TransBioLine project (No 821283).

Project. The sum of all activities carried out in the framework of the Grant Agreement.

Consortium. The TransBioLine Consortium, comprising the above-mentioned legal entities.

Consortium Agreement. Agreement concluded amongst TransBioLine participants for the implementation of the Grant Agreement. Such an agreement shall not affect the parties' obligations to the Community and/or to one another arising from the Grant Agreement.

Publishable Summary

The Translational Safety Biomarker Pipeline (TransBioLine) is an international project with 27 partners in 10 countries. The Central Biomaterial Bank Charité (ZeBanC) is responsible for the biobanking including coordination of sample-related processes, the provision of sample kits, storage and distribution of samples to the analysis laboratories. In addition, ZeBanC contributes significantly to the harmonization of data management and data flow.

Projects of this dimension suffer from a considerable heterogeneity of partners and contributors. As every recruitment site and analysis partner has their own standard procedures and workflows it is very important to harmonize and standardize the fundamental processes.

The first major deliverable for work package 7 is to harmonize standard operating procedures (SOPs) for sample collection, processing, shipment, aliquoting and shipment. This harmonization should lead to samples and related data of high quality and comparability.

1. Introduction

The aim of this project deliverable is to harmonize standard operating procedures (SOPs) for sampling, processing, interim storage and shipment. This step is of high importance due to the heterogeneity of the 27 project partners located in 10 different countries and to the variety of biomarkers to be analysed. Differences in the pre-analytical procedures can affect sample quality and might thus cause analytical data reflecting the source of the samples rather than biological variations. The standardization of the mentioned processes is intended to guarantee the consistent sample and data quality for biomarker analysis.

2. Methods

2.1 Retrieval of information

2.1.1 Request of information for analyses

In order to assess the requirements of the different TBL collections, a survey for collection of process and sample details was initiated. All information on the collection strategy (sample type, number) and sample handling such as centrifugation parameters etc. were compiled. In this context, questions were also asked about the biomarkers to be analyzed. For a better understanding, the biomarker companies were also involved in order to include their requirements (sample quantities and quality parameters, if applicable) into the survey.

2.1.2 Information study protocols/SOPs

Organ-related work packages (WP1 to WP 5) and EFPIA partners were asked to contribute with their own protocols/SOPs.

2.2. Harmonization of SOPs

2.2.1 SOP for sample collection and processing

The ZeBanC elaborated common procedures and adapted similar workflows where possible. However, minor differences due to local settings and standards will not be avoidable.

After retrieving all information the ZeBanC created a harmonized protocol focusing primarily on the requirements of the partners for biomarker analyses to cover optimal sample conditions. The local staff will be trained according to this protocol in order to allow compliance to the agreed protocol.

3. Results

3.1 SOP for sampling, processing and intermediate storage

The finalised document KP-08_04.29a SD TransBioLine Processing Manual version 1 was uploaded to the TransBioLine SharePoint for availability of all TBL partners. Additionally, the document was send on request to the staff member without direct access to SharePoint. To facilitate implementation especially in daily practice, the ZeBanC produced movies with detailed explanation of all the important processes such as sampling, processing of samples and shipment and uploaded the movies to SharePoint. If desired, the ZeBanC will organise a video conference training the staff. This way, the staff members can follow a specific workflow to retrieve all required information.

Deviations from the protocol such as different centrifugation conditions, varying temperatures have to be documented in the eCRF and/or on the sample data sheet.

To allow easier identification of the different sample types at the recruitment sites, a color-coding system for the tubes is used. Different colours of the tube lid for different types of materials facilitate easier and safe aliquoting (e.g. orange for serum, green for EDTA-plasma, blue for the 2x centrifuged plasma (defined as liquid biopsy) and yellow for urine). The volume and number of aliquots to be collected need to be confirmed by each work package. For serum and urine measured at MLM 5 ml barcoded tubes are used which will be labelled accordingly.

Each sample kit that is provided is accompanied with sample data sheets. This document carries the respective barcodes and sample IDs and the general sample information should be collected e.g. TBL patient ID (provided by eCRF or ABX-CRO), date and time of sampling etc. This data sheet should assist in collecting the necessary data, which then can be entered into the respective eCRF.

3.2. SOP for storage

The ZeBanC already has standard procedures for sample storage. All samples from TransBioLine (TBL) will be stored at -80°C and documented in our laboratory information and management system (LIMS, here CentraXX). 2D barcoded tubes (containing 2D barcode, 1D barcode and human readable number) will be used for easier and unique identification. Samples will be stored in defined spaces, are tracked and easily detectable from freezer to rack and single tube. A standard procedure for the removal/destruction of samples and data in case of consent withdrawal has still to be defined. For this, all the informed consents from the organ work packages have to be available and the responsibilities and contacts have to be defined.

3.2.3 SOP for shipment

After sample collection, each organ work package and/or location, which is not using the eCRF has to send a data export to ABX-CRO. There the data will be checked and a TBL-patient-ID will be allocated to each patient. This ID is very important for patient

and sample identification within the TransBioLine project. For work packages using the eCRF, the patient ID will be generated automatically.

Prior to each sample shipment to the biobank, each partner has to upload the data of the samples to their specific ownCloud folder (contact ITTM) and inform the biobank and ABX-CRO. ABX-CRO will cleanse the data and clear for usage. The biobank will check the data and confirm shipment. Only after this confirmation, the samples can be transferred to the biobank.

The SOP for shipment is included into the document KP-08_04.29a SD TransBioLine Processing Manual version 1. Additionally, a video on packaging for shipment was uploaded to SharePoint. It describes the general handling and packaging of biological samples, which are assigned to UN3373 regulations. As every site is using their own transport service small deviations may apply.

4. Discussion

(Include if relevant deviations from DoA and contingency plans)

Biobanking in a heterogeneous research consortium requires harmonization in order to generate reliable and reproducible data. Different data structures, different local infrastructures, different qualification of staff members and diverse procedures must be brought to a comparable level.

The harmonization of SOPs depends not only on the different organ-specific work packages and their recruitment sites, but also on the respective partner for biomarker analyses. The harmonized SOP is based on a compilation of the various study protocols and the requirements of the partners for biomarker analyses. Deviations from these protocols have to be documented. Minor changes such as the use of different tubes due to their own vendors were tested and implemented.

5. Conclusion


The respective documents and videos for sample collection, processing, storage and shipment have been created and accessible for everyone within the TransBioLine consortium. If desired, the ZeBanC will organise a video conference training the staff. Whenever questions occur, the ZeBanC staff is happy to support.

Annexes

Annex 1

KP-08_04.29a SD Transbioline Processing Manual

KP-08_04.30 VA TBL-Prozesse intern

TransBioLine: Processing Manual Harmonized SOPs

Contents

- (3) Request of sample kits
- (4) Sample collection containers
- (5) Sample kits - I
- (6) Sample kits - II
- (7) Sample data sheet – general information
- (8) Sample data sheet – patient information
- (9) Sample data sheet – sample information
- (10) Sample data sheet - barcodes
- (11) Sample processing
- (12) Serum (SR)
- (13) EDTA-plasma (EP)
- (14) Liquid biopsy (LB)
- (15) EDTA-Plasma/Liquid biopsy
- (16) Urine (UR)
- (17) Shipping
- (18) Packing UN3373
- (19) Contact/ Shipping address

Request of sample kits

3

Order form TransBioLine sample kits

Please fill in the form and send via email to zebanc-labor@charite.de

Please send us the following sample kits:

1)

Work package	<input type="checkbox"/> WP1 – DIKI	<input type="checkbox"/> WP3 – DIP1	<input type="checkbox"/> Healthy volunteers (HV)
	<input type="checkbox"/> WP2 – DILI	<input type="checkbox"/> WP4 – DVI	
Location	<input type="checkbox"/> CDG	<input type="checkbox"/> MUC	<input type="checkbox"/> SLM
	<input type="checkbox"/> TXL	<input type="checkbox"/> KEF	<input type="checkbox"/> LID
	<input type="checkbox"/> EFP	<input type="checkbox"/> LPL	<input type="checkbox"/> NCL
	<input type="checkbox"/> PST	<input type="checkbox"/> NQT	<input type="checkbox"/> ZRH
	<input type="checkbox"/> BCN	<input type="checkbox"/> AGP	<input type="checkbox"/>

2)

Number of kits
(patients)

3)

Contact person for further enquiry: _____
Telephone number/ e-mail: _____
Date of order : _____
Signature: _____
.....

Delivery sent (via courier service)

5)

Date : _____
ZeBanC signature: _____
.....

Receipt (please e-mail document to zebanc-labor@charite.de)
Date of receipt sample kits: _____
Signature : _____
Note : _____

Please request sample kits via the provided order form and send to zebanc-labor@charite.de

- 1) **Work package** - please, tick respective workpackage
- 2) **Location** - please, tick location code, list will be provided
- 3) **Number of kits** – number of racks and sets (patients) in brackets (will differ between work packages)
- 4) **Contact** – please, fill in contact details in case of questions
- 5) **Receipt** – The request sheet will be enclosed in the package. Please fill in date of sample reception.


Note: Please order at least 1 week in advance!

Sample collection containers

- Master samples

Plasma

- BD Vacutainer™ Hemogard Closure Plastic K2-Edta Tube, 10 ml (Brand: BD 367525)
- VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 lavender cap-black ring, PREMIUM (Brand: Greiner 456036)

Serum

- BD Vacutainer® Plastic Serum Tube with Red BD Hemogard™ Closure, 10ml (Brand: BD 367896)
- VACUETTE® TUBE 9 ml CAT Serum Separator Clot Activator, gold cap-gold ring (Brand: Greiner 455034)


Urine

- BD Vacutainer™ Plastic Urinalysis Tubes (Brand: BD 364915)

If other, please specify on data sheet

Sample kits - I

Bag with 5 ml barcoded tubes (MLM)

Box for 5 ml barcoded tubes

Sample data sheet

Rack with 0.8 ml barcoded tubes

- One sample kit usually contains tubes for several samplings (several sets)

Example WP2

- 1 Rack 0.8 ml barcoded tubes (8 patients, one patient per row)
- 1 bag with 5 ml barcoded tubes (8, one per patient)
- Several empty racks (0.8 ml) and boxes (5ml) for storage in freezer

Sample kits - II

6

- The ZeBanC will provide sample kits with 2D barcoded tubes and sample data sheets for each WP. The number of tubes will differ, but the layout will be the same:

- Serum**

0.8 ml tubes

5 ml tubes

- EDTA-Plasma**
(1x centrifuged)


- Liquid Biopsy**
(2x centrifuged, cell free plasma)

- Urine**

Infobox Barcodes

5 ml

0.8 ml

Sample data sheet – general information

1)

Sample Kit-ID	<<<(1)>>>	<<<(d1)>>>	
Rack ID	<<<(r1)>>>	<<<(s1)>>>	
TBL Patient ID*	_____		
Pro Euro DILI ID	_____		
Site name**	<input type="checkbox"/> NQT <input type="checkbox"/> AGP <input type="checkbox"/> BRN <input type="checkbox"/> LIS <input type="checkbox"/>		
Patient population**	_____		
Visit	<input type="checkbox"/> V1 <input type="checkbox"/> V1b <input type="checkbox"/> V1c <input type="checkbox"/> V2 <input type="checkbox"/> V3 <input type="checkbox"/> V4 <input type="checkbox"/> V5		
Sample type	Serum [SR]	EDTA-Plasma [EP]	Liquid Biopsy [LB]
Sampling date	20 _ _ / _ _ / _ _		
Sampling time	_ _ : _ _	_ _ : _ _	_ _ : _ _
Master sample ID*	_____		
Volume master sample	____ ml	____ ml	____ ml
Start 1. Centrifugation	_ _ : _ _	_ _ : _ _	_ _ : _ _
Time of freezing	_ _ : _ _	_ _ : _ _	_ _ : _ _
Storage location 0.8 ml tubes	Rack ID: _____ / Freezer: _____		
Storage location 5.0 ml tubes	Rack ID: _____ / Freezer: _____		
Sampling by***	Processing by***		


* Generated in the eCRF mask or will be provided by the Biobank, ** Code can be found in the document "WP2 recruitment site and patient population matrix", *** only for site-specific questions, will not documented in the eCRF nor in the Biobank

1) General information for Biobank

2) Processing

- a) Overview processing
- b) Aliquoting scheme

- Each sample kit will include tubes for several patients/visits (here one patient per row = row A comprises 4 orange, 6 green and 2 blue tubes, next row, row B =next patient)
- Please check, which row(s) you will need for your patient!

Sample data sheet – patient information

Please fill in the following information for each patient/visit:

1) TBL Patient ID*	_____		
2) Pro Euro DILI ID	_____		
3) Site name**	<input type="checkbox"/> NQT <input type="checkbox"/> AGP <input type="checkbox"/> BRN <input type="checkbox"/> LIS <input type="checkbox"/> _____		
4) Patient population**	_____ _____		
5) Visit	<input type="checkbox"/> V1 <input type="checkbox"/> V1b <input type="checkbox"/> V1c <input type="checkbox"/> V2 <input type="checkbox"/> V3 <input type="checkbox"/> V4 <input type="checkbox"/> V5		
Sample type	Serum [SR]	EDTA-Plasma [EP]	Liquid Biopsy [LB]
Sampling date	20 ____ / ____ / ____		
Sampling time	____ : ____	____ : ____	____ : ____
Master sample ID*	_____		
Volume master sample	____ ml	____ ml	____ ml
Start 1. Centrifugation	____ : ____	____ : ____	____ : ____
Time of freezing	____ : ____	____ : ____	____ : ____
Storage location 0.8 ml tubes	Rack ID: _____ / Freezer: _____		
Storage location 5.0 ml tubes	Rack ID: _____ / Freezer: _____		
Sampling by***	Processing by***		

* Generated in the eCRF mask or will be provided by the Biobank, ** Code can be found in the document " WP2 recruitment site and patient population matrix". *** only for site-specific questions, will not documented in the eCRF nor in the Biobank

- 1) **TBL (Transbioline) Patient ID** – will be created by eCRF or if data transfer by ABX-CRO
- 2) Study specific ID – only if applicable e.g. DILI
- 3) **Site name** - specific TBL site code, a list for each WP will be provided
- 4) **Patient population**- specific population code for each WP, a list for each WP will be provided
- 5) **Visit** – Visit number

Sample data sheet – sample information

Please fill in the following information for each patient/visit:

TBL Patient ID*	_____		
Pro Euro DILI ID	_____		
Site name**	<input type="checkbox"/> NQT <input type="checkbox"/> AGP <input type="checkbox"/> BRN <input type="checkbox"/> LIS <input type="checkbox"/> _____		
Patient population**	_____ _____ _____		
Visit	<input type="checkbox"/> V1 <input type="checkbox"/> V1b <input type="checkbox"/> V1c <input type="checkbox"/> V2 <input type="checkbox"/> V3 <input type="checkbox"/> V4 <input type="checkbox"/> V5		
Sample type	Serum [SR]	EDTA-Plasma [EP]	Liquid Biopsy [LB]
Sampling date	20 ____ / ____ / ____		
Sampling time	____ : ____	____ : ____	____ : ____
Master sample ID*	_____		
Volume master sample	____ ml	____ ml	____ ml
Start 1. Centrifugation	____ : ____	____ : ____	____ : ____
Time of freezing	____ : ____	____ : ____	____ : ____
Storage location 0.8 ml tubes	Rack ID: _____ / Freezer: _____		
Storage location 5.0 ml tubes	Rack ID: _____ / Freezer: _____		
Sampling by***	Processing by***		

* Generated in the eCRF mask or will be provided by the Biobank, ** Code can be found in the document " WP2 recruitment site and patient population matrix". *** only for site-specific questions, will not documented in the eCRF nor in the Biobank

- 6) **Sampling Date** – format yyyy/mm/dd
- 7) **Sampling time** – time point when sample was taken
- 8) **Master sample ID** – automatically generated by eCRF mask or provided by the Biobank
- 9) **Volume master sample** – volume in millilitre (ml) in original sampling device before centrifugation
- 10) **Start 1. Centrifugation** – time point 1st centrifugation (important for data about sample quality)
- 11) **Time of Freezing** – time point of freezing at -80°C
- 12) **Storage location** – temporary freezer location of rack at recruitment site
- 13) **Sampling by/ Processing by** – name of person, who sampled the biospecimens and name of person who performed centrifugation, aliquoting

Sample data sheet - barcodes

10

96 well ROW A / SET 1 0.8 ml 2D Barcode Tube			Single Tube No. WP2_S1 5 ml 2D Barcode Tube		
Serum [SR] 0.5 ml			Serum [SR] 0.5 ml		
SR_1	FR18305527		SR_5	B001666455	
SR_2	FR18305549				
SR_3	FR18305534				
SR_4	FR18305542				
Plasma [EP] 0.5 ml					
EP_1	FR16971733				
EP_2	FR16971729				
EP_3	FR16971731				
EP_4	FR16971728				
EP_5	FR16971792				
EP_6	FR16971706				
Liquid Biopsie [LB] 0.5 ml					
LB_1	FR22461488				
LB_2	FR18305528				


Please cross out unused tubes.
Collect unused tubes and send them back to the biobank.

- Overview about individual tubes with human readable number and barcode
- Please cross out all the tubes you don't use, collect them and send them back to biobank.
- Usage of 5 ml tubes depend on WP, but are labeled according to sample data sheet (here SR_5)
- Always compare numbers with barcodes on tubes!

Sample processing

11

- Invert tubes according to manufacturer's instruction

Serum (SR)

12

- (1) Collect blood for serum
- (2) Invert tube **4-6 times**
- (3) Allow to coagulate for 30-60 min at room temperature
- (4) Centrifuge (within 2 hours) at **2000xg, 10 min, 4-8 °C (or RT)**
- (5) Carefully aliquot 0.5 ml supernatant (serum) into each 2D barcoded tube (**0.8ml orange cap, 5ml white cap**)
- (6) Freeze immediately and store at -80°C

supernatant (serum)

clot

- (1) Collect whole blood in EDTA-tubes
- (2) Invert tube 8-10 times
- (3) Store on wet ice in upright position
- (4) Centrifuge (within 2 hours) at **2.000xg, 10 min, 4-8 °C**
 - a) Aliquot 0.5 ml into 2D barcoded tubes (**EP, 0.8ml green cap**)
- (5) Freeze immediately and store at -80°C

- (1) Collect whole blood in EDTA-tubes
- (2) Invert tube 8-10 times
- (3) Store on wet ice in upright position
- (4) Centrifuge (within 2 hours) at **2.000xg, 10 min, 4-8 °C**
- (5) Carefully collect supernatant without disturbing the cell pellet and transfer to **new tube**
- (6) Centrifuge **new tube** at **10.000xg, 10 min, 4-8 °C**
 - a) Aliquot 0.5 ml into 2D barcoded tubes (**LB, 0.8 ml blue cap**)
- (7) Freeze immediately and store at -80°C

Only if using one master sample (collection tube) for EDTA-plasma and LB, you should continue here:

- (1) Collect whole blood in EDTA-tubes
 - (2) Invert tube 8-10 times
 - (3) Store on wet ice in upright position
 - (4) Centrifuge (within 2 hours) at **2.000xg, 10 min, 4-8 °C**
 - (5) Carefully collect **1.2 ml** of supernatant without disturbing the cell pellet and transfer to **new tube**
 - (6) Aliquot 0.5 ml into 2D barcoded tubes (**EP, 0.8 ml green cap**)
-
- (7) Centrifuge **new tube** at **10.000xg, 10 min, 4-8 °C**
 - a) Aliquot 0.5 ml into 2D barcoded tubes (**LB, 0.8 ml blue cap**)
 - (8) Freeze immediately and store at -80°C

- (1) Collect urine
- (2) Allow to cool down
- (3) Aliquot 0.5 ml into 2D barcoded tubes (**UR,0.8 ml yellow cap**)
- (4) Freeze immediately and store at -80°C

- The ZeBanC in Berlin has to be informed about the date of sample shipment in advance.
- The Biobank has to confirm its acceptance to receive a shipment of samples before the sample transfer.
- Samples have to be transferred frozen to the ZeBanC
- For transfer of frozen samples to the ZeBanC in Berlin make sure that the samples will be under dry ice all the time.
- For shipment use a carrier or shipping agent that guarantees delivery of samples in frozen status to the ZeBanC
- **Please note that delivery of samples to the ZeBanc should be from Monday to Thursday (9:00 to 16:00).**

Blood and its components have to be packed according to **UN3373** – biological substance category B

1. Wrap rack with aliquots (2D barcoded tubes) in fleece (A)
2. Put this into bag (B) and seal
3. Put racks on dry ice in styrofoam box (C)
4. Label box with address (always include name and phone number of sender and recipient) and UN 3373 sticker (D)
5. **Attention:** if air cargo, additional label for dry ice (E) and arrows (F) are needed

Usually, all provided by courier service!

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